





TACTICAL PACRAC™ (TPR)

The Tactical PacRac™ (TacRac) is a lightweight equipment platform compatible with most NATO litters (STANAG 2040) and can be configured to secure a wide array of medical equipment.

Art. no	274-M
Width	48 cm
Folded Height	7,6 cm
Deployed Height	33 cm
Depth	34,3 cm
Weight	4,5 kg (unequipped)
NSN	6530-66-166-0063



Ideal for field hospitals, the TacRac serves as a complementary equipment table to our Compact Treatment Station. Designed for mobility, the TacRac can be rapidly deployed or removed as needed. Its folding legs allow for easy stowage when not in use, while the folding sides enable packing into hard cases or load-out bags.

The TacRac can be configured to secure various medical devices and follow the litter through the MedEvac chain. It provides a lightweight, modular equipment table for field treatment, keeping medical devices organized and secure during treatment and CasEvac. Constructed from 6061-T6 aluminum, the TacRac is built to withstand harsh environments.