



## Defibrillator - Mount

### AED Mount Wall iNTRAXX

The AED mount is for iNTRAXX wall track mounting. Can mount into vertical or horizontal track. The AED can be locked into this mount with or without its protective carry case.

	IMPERIAL	METRIC
<b>Part Number</b>	F90241	
<b>Dimensions</b>	W 10.2 in D 9 in L 14.4 in	W 260 mm D 229 mm L 367 mm
<b>Weight</b>	6.9 lb	3.1 kg
<b>Load Capacity</b>	7.7 lb	3.5 kg
<b>Crash Tested</b>	10g, 20g, 26g	

- Supplied complete with bolts and spacers to suit each AED
- Can be retained in the mount with or without the AED Case always refer to instructions
- iNTRAXX wall mount only
- Refer to separate customer drawings for mounting details
- Defibrillator not included
- Please refer to separate operators manual
- Ferno is not responsible of testing claims of OEM mounts
- Currently known devices, Zoll AED plus & Mindray BeneHeart D1

#### KEY

- (W) Width is measured as 'Y' axis in relation to mounting to the track
- (D) Depth is measured as 'Z' axis in relation to mounting to the track
- (L) Length is measured as 'X' axis in relation to mounting to the track

