





CUSTOMER DRAWING

MOUNT IN TRACK


MINDRAY BENEHEART
VERSION D1

MINDRAY BENEHEART
VERSION C1A



NOTES:

- APPROX. WEIGHT:1,8 KG

GENERAL TOLERANCES COMPLY WITH STN ISO 2768-mK		THIS DRAWING IS THE PROPERTY OF FERNO SLOVAKIA S.R.O.. IT IS NOT TO BE COPIED OR DUPLICATED IN ANY MANNER, AND SHALL NOT BE SUBMITTED TO OUTSIDE PARTIES FOR EXAMINATION WITHOUT OUR KNOWLEDGE AND CONSENT. IT SHALL BE USED AS A MEANS OF REFERENCE TO WORK FURNISHED BY THIS COMPANY ONLY.		
DRAWN BY A.TURIAK	DATE 27.05.2026	 FERNO-SLOVAKIA, s.r.o. Bosaca 893, 913 07 Slovak Republic		
MATERIAL		FIRST ANGLE PROJECTION 	SCALE: 1:5	SHEET: 2 OF 2
FINISH		DESCRIPTION MINDRAY MOUNT - TRACK MOUNTED		
REFERENCE NONE		SIZE A3	PART NO. 048-0268	REVISION 0A