


# CUSTOMER DRAWING


MINDRAY BENEHEART  
VERSION D1

MINDRAY BENEHEART  
VERSION C1A

4xM6x16(DIN7991)  
(included)




## MOUNT



### NOTES:

- APPROX. WEIGHT:1,8 KG

|                                                          |                           |                                                                                                                                                                                                                                                                                              |                             |                         |
|----------------------------------------------------------|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|-------------------------|
| GENERAL TOLERANCES COMPLY WITH<br><b>STN ISO 2768-mK</b> |                           | THIS DRAWING IS THE PROPERTY OF FERNO SLOVAKIA S.R.O.. IT IS NOT TO BE COPIED OR DUPLICATED IN ANY MANNER, AND SHALL NOT BE SUBMITTED TO OUTSIDE PARTIES FOR EXAMINATION WITHOUT OUR KNOWLEDGE AND CONSENT. IT SHALL BE USED AS A MEANS OF REFERENCE TO WORK FURNISHED BY THIS COMPANY ONLY. |                             |                         |
| DRAWN BY<br><b>A.TURIAK</b>                              | DATE<br><b>27.05.2026</b> |  <b>FERNO-SLOVAKIA, s.r.o.</b><br>Bosaca 893, 913 07 Slovak Republic                                                                                                                                    |                             |                         |
| MATERIAL                                                 |                           | FIRST ANGLE PROJECTION<br>                                                                                                                                                                              | SCALE:<br><b>1:5</b>        | SHEET:<br><b>2 OF 2</b> |
| FINISH                                                   |                           | DESCRIPTION<br><b>MINDRAY MOUNT - WALL MOUNTED</b>                                                                                                                                                                                                                                           |                             |                         |
| REFERENCE<br><b>NONE</b>                                 |                           | SIZE<br><b>A3</b>                                                                                                                                                                                                                                                                            | PART NO.<br><b>048-0269</b> | REVISION<br><b>0A</b>   |