



XTRIC8 L

VERSATILITY MEETS URGENCY.


The Ferno XTRIC8 Rescue Sled is a versatile and robust extrication solution, engineered for speed, reliability, and adaptability across a wide range of operational environments. Whether operating in lowland terrain, mountainous regions, confined spaces, or conducting high-angle extractions, the XTRIC8 delivers consistent, high-performance results when every second is critical.

Compact and durable, the XTRIC8 is purpose-built for rapid deployment in the most demanding conditions. Its intelligent design streamlines casualty packaging and evacuation, while maintaining the tactical flexibility and traditional functionality required by military and defence personnel.

Key Features

- » **Pre-rigged for rapid deployment** and immediate use straight from the carry bag. No assembly required.
- » An **integrated foot restraint** ensures maximum casualty security during vertical lifts and high-angle extractions.
- » Colour-coded, **four-point quick-connect lifting bridle** is pre-set for fast and efficient rigging.
- » Straps and buckles are designed for **ease of use with gloved hands**, ensuring swift and secure fastening and unfastening.
- » **Heavy-duty carry bag** features four large compartments for additional gear, and an adjustable backpack harness for easy transport.
- » XTRIC8's user-friendly design promotes **intuitive operation**, **reduces training time** and **enhances operational readiness** in high-pressure scenarios.





SPECIFICATIONS

XTRIC8 L (TACTICAL VERSION)

NATO NSN	6530-66-168-8188
Manufacturer's Product Code	FWE XTRIC8 L
Inclusions	XTRIC8 L Rescue Sled pre-rigged with internal patient foot strap, external restraints, external straps and handles
	1 x Maillon Rapide for vertical lifting
	1 x Colour coded 4-Point Lifting Bridle
	1 x Storage Bag with backpack straps
Length	2540 mm
Width - Widest Point	914 mm
Weight - XTRIC8 L only (rigged)	6.5 kg
Weight - XTRIC8 L (rigged) and Bag	8.5 kg
Load Limit - XTRIC8	350 kg - HORIZONTAL / 250 kg - VERTICAL
XTRIC8 LIFTING BRIDLE	
Weight - Lifting Bridle	0.9 kg
Load Limit - Lifting Bridle	500 kg
Bridle Karabiners	VYO N-2853G - Aluminium Alloy (SAE 7075) Triple Action Karabiner
Bridle Karabiner Dimensions	124 mm x 80 mm / Gate Opening: 22 mm
Bridle Karabiner Ratings	Minimum Tensile Strength: 23 kN / Minimum Gate Strength: 16 kN
Bridle Karabiner Standards	ANSI Z359.12-2019; CSA Z259.12-16; EN 362:2004

- Ferno and manufacturers of components reserve the right to change specifications without notice.
 - General specifications are rounded to the nearest whole number. Metric conversions are calculated before rounding the English measurements.
- For more information contact Ferno Customer Service or your local distributor.