



TRANSFER MATTRESS

STRETCHER TRANSFER MATTRESS COMPATIBLE WITH ALL NATO STRETCHERS WITH 4-POINT RESTRAINTS

Specifically designed for use with NATO STANAG 2040 stretchers (with handles folded in) and for CasEvac operations. It can be easily folded or rolled up for convenient storage. The underside of the mattress features a glass fiber surface that reduces friction, enabling smooth transfers and easy sliding of both the mattress and the patient from bed to stretcher. It is also suitable for use in confined spaces where a stretcher may not fit.

The mattress handles support loads up to 250 kg, with fully welded seams for maximum hygiene and easy cleaning. The water-resistant padding adds durability. Additionally, it can function as a standalone stretcher for short CasEvac transport. When paired with the Step by Step system, this solution provides exceptional ergonomics and load-bearing capacity, supporting up to 250 kg. The mattress is approved according to European Standard EN 1865.

Width	65 cm (47 folded)
Height	3,5 cm
Lenght	188 cm
Weight	4,5 kg
Art. Nr	GE21271008881
Supplier Nr	271008880
Certifications	MDR 2017/745 Class 1



KEY FEATURES

- MDR 2017/745 Class 1 concerning Medical Devices.
- European Standard EN 1865.
- Heat/Cold resistance: +70c / -30c
- Fits all NATO STANAG 2040 Stretchers/Litters
- Low IR reflection, tested according to FSD 7403 with an average value of 30-45 nm