



# THERMO SAVER HYPOTHERMIA BAG

*Thermosaver is a rescue bag that is small, lightweight, quick to use, and effective for the evacuation of hypothermic and injured persons. It has been developed in collaboration with Norwegian experts in rescue operations and emergency medicine, as well as the Norwegian Defence Research Establishment.*

<b>Art. Nr</b>	MQ20016045
<b>Certifications</b>	EN-1865
<b>Weight</b>	2600 g
<b>NSN</b>	6532.2516.35297



An “All-in-one” hypothermia rescue bag that replaces windproof cloths, wool blankets, bubble wrap, and carrying slings. The rescue bag has a water-resistant outer fabric, a polyester-based insulation material in the middle, held in place by a thin polyester sheet. In the lying area and under the head, the insulation is double thick. A water-absorbing moisture barrier (ASAP/Gemini) should be used closest to the patient.

The burrito principle is used to wrap the patient, meaning the rescue bag is flexible and fits all body sizes. The bag, with the inserted moisture barrier, is rolled up to enable the “log-roll” method for moving the patient onto the rescue bag. Once the folding is completed around the patient, the rescue bag is secured using Velcro fasteners.

Tested to EN-1865 standards with a lifting strength of at least 250 kg for each carrying handle.



## RELATED PARTS

Part NO	Description	Price	UOM	Quantity	Unit
20016044	First Aid Thermosaver Water Repellent (Yellow)	9799	pcs	1	pcs
20016045	First Aid Thermosaver M Water Repellent (Grey)	9799	pcs	1	pcs
20016046	First Aid Thermosaver MJK Water Repellent (Grey)	9363	pcs	1	pcs
20016047	First Aid Thermosaver Hood Water Repellent (Yellow)	1842	pcs	1	pcs
20016048	First Aid Thermosaver Hood M Water Repellent (Black)	1842	pcs	1	pcs
20016049	First Aid ASAP Jona for Thermosaver	2169	pkg	1	pcs
20016050	First Aid ASAP Gemini for Thermosaver	1733	pkg	1	pcs
20016051	First Aid Thermosaver Compression Bag (Yellow)	1797	pcs	1	pcs
20016052	First Aid Thermosaver M Compression Bag (Grey)	1797	pcs	1	pcs
20006411	Ready Heat Active Warming Package 1-panel 41x46 cm	259	pcs	1	pcs
10024593	Ready Heat Active Warming Package 6-panel 86x61 cm	5668	pkg	12	pcs