

Product Datasheet - Thermosaver Compression bag

Manufacturer: Thermosaver AS

Product Description

Compression bag made for Thermosaver Hypothermic Rescue Bag.

- The compression bag has the same material as the outer layer of the rescue bag itself, water-repellent polyester and comes in both hi-vis yellow and grey.
- The bag has three tightening straps to compress the bag to the maximum.
- There are handles at both ends to get the rescue bag out and in use as quickly as possible.
- The compression bag also has a shoulder strap that makes it possible to have the rescue bag on your back.
- The bag can be used several times and can be recycled.



NATO Stock Number (NSN): 6545251635298

© Thermosaver 2023

Page 1 of 2

Document Nr. and Revision: PS-032 Rev. 4

Release date: 31.08.2023

Production Specification

Compression bag	
In Use L/W/H cm	70/49/1
Weight g	300
Colour	Hi-vis Yellow / Grey
Fabric	100 % polyester
Outer skin fabric watercolumn	8000 mm
Shrinkage wash 60 deg. C	<3 %
Biocompatibility lining	Oeko Tex 100
Visibility	ISO 20471