





PELICAN GOLDEN HOUR MEDIC CONTAINER

The Golden Hour™ Medic Container is a nimble and tactical option for delivery of whole blood to the point of injury within a crucial window of time. Capable of keeping two liters of whole blood cool for 24 hours under extreme conditions and up to three days under optimal conditions.

Color	MultiCam
Width	35.56 cm (14 in.)
Height	15.24 cm (6 in.)
Depth	22.86 cm (9 in.)
Weight	2.722 kg (6 lbs)



The MOLLE attachment points allow the reusable carrier to be mounted on medic's kit and reduced footprint increases mobility over the traditional cumbersome containers.

Utilizing Thermal Isolation Chamber System (TIC™) technology and vacuum insulation panels, the Series-S Medic Container will keep the payload isolated from exposure to extreme temperatures with a tight, protective thermal seal for a pre-determined temperature range and duration.

Safely deliver whole blood and blood supplies from anywhere in the world — to where they need it most. The containers are easy to assemble, reusable parcels that are constructed of reusable non-toxic materials.



SPECS

- Stores two liters of whole blood
- Keep blood cool for 24 hours under extreme conditions, up to 3 days under optimal conditions
- Designed and tested for harsh environments
- Reusable carrier mounts on medic's kit
- GMC pouch includes accessory storage options