



PELICAN GOLDEN HOUR ONE-V CONTAINER

The Golden Hour™ ONE-V Container is an ultra-lite blood transport capable of carrying one half liter of whole blood and keeping the stored blood cool for up to 48 hours, regardless of weather conditions.

| | |
|---------------|------------------|
| Color | MultiCam |
| Width | 15.24 cm (6 in.) |
| Height | 20.32 cm (8 in.) |
| Depth | 12.7 cm (5 in.) |
| Weight | 0.907 kg (2 lbs) |



This container easily transports a single unit of whole blood to any location. This extremely light-weight container will securely and reliably transport blood and platelets to wherever it is needed most.

Utilizing Thermal Isolation Chamber System (TIC™) technology and vacuum insulation panels, the Golden Hour™ ONE-V Container will keep the payload isolated from exposure to extreme temperatures with a tight, protective thermal seal for a pre-determined temperature range and duration.

Safely deliver whole blood and blood supplies from anywhere in the world — to where they need it most. The containers are easy to assemble, reusable parcels that are constructed of reusable non-toxic materials.



SPECS

- Ultra-light storage and transport of one half liter of whole blood
- Keep blood cool for up to 48 hours, regardless of weather conditions
- Dismount for up to 12 hours without external power
- Equip to make blood in far forward areas accessible
- Reusable, rugged design for easy cleaning and replacement